	SEROCON R&D BIOTECHNOLOGY HEALTH CHEMICAL INDUSTRY AND TRADE INC	Document Code	PRP.12
		Effective Date	01.11.2022
		Rev. No / Date	02/22.12.2025
		Page	1/2
Therapeutic Drug Monitoring External Quality Control Prospectus			

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LIMITATONS

1. This product must not be used after the expiration date.
2. A syringe must not be used for reconstituting the sample.
3. If there is any indication of microbial contamination or excessive turbidity in the reconstituted product, the vial must be discarded.

TEST TIMES

According to the program code, the box contains 4-month (TDM4) samples. Each sample is labeled with the cycle month for entry of results into the SEROCON Portal via www.serocon.com. Samples should be analyzed on the dates specified for the corresponding cycle month according to the program.

REPORTING OF RESULTS

Results must be entered into the SEROCON Portal via www.serocon.com using a username and password no later than the last day of the relevant cycle month. When entering results, each institution uses the program opened for its own device and method. If a device change occurs during the program, SEROCON must be informed so that necessary adjustments can be made. In case of a device change, participants leave the data from the old device as of the date of change and enter new results according to the defined new program schedule. Results are released starting from the second week of the month following the cycle. The system does not allow entry of results outside the limit values for each parameter. Results that should be entered with ">" or "<" must be entered without these symbols. Participants can view their reports via the SEROCON PORTAL. The program is based on the periodic presentation of results obtained for each parameter from the analysis of samples with unknown concentrations. For each parameter, all participants' results are calculated according to ISO 13528 requirements, and Z or Z' scores (SDI) are generated. Values entered by other participants can be viewed on graphs without revealing participant names, in accordance with confidentiality principles. To help participants track their performance, values from previous cycles are provided as Levey-Jennings charts for each parameter at the sample analysis periods.

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Automatic pipette
Pipette tip
Distilled or deionized water

ABBREVIATIONS

PC	:	Program code
CONT	:	Contents
	:	Temperature Limitation
LOT	:	Lot Number
IVD	:	In Vitro Diagnostic Medical Device
CE	:	European Conformity

Note: You need to obtain the current documents by logging into SEROCON PORTAL via www.serocon.com website. Shared Documents;

- ✓ SEROCON DKD programs protocol
- ✓ EQC Programs Work Schedule
- ✓ Program Prospectuses

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SAMPLE WORK SCHEDULE

Working Months**	Sample Lots	Proposed Study Date	Results Entry Deadline
February	1702602	21.02.2026	28.02.2026
May	1702605	23.05.2026	31.05.2026
August	1702608	22.08.2026	31.08.2026
November	1702611	24.11.2026	30.11.2026

**The sample vials display the Program Code and LOT number. The last two digits (XX) of the LOT number indicate the sample number for the relevant cycle month. The user performs the analysis by referring to the XX number for the corresponding cycle month.

CONTACT US

SEROCON R&D BIOTECHNOLOGY HEALTH CHEMICAL INDUSTRY AND TRADE INC

Fevzi Çakmak District, 10739. Street No: 16-42050, Karatay, Konya, Türkiye

Tel: 0850 303 6644

Fax: 0332 353 8488

info@serocon.com

www.serocon.com